


SHOP MICROSCOPES



For printing, machining, inspection, measurements

Model No.	M1752	M1754
Magnification(Field of View)	40x(4.5mm)	100x(1.8mm)
Objective Lens	4x	10x
Focusable Eyepiece	10x with scale	10x with scale
Measuring Range	3.0mm / 0.1"	1.0mm / 0.05"
Working Distance	15mm	5.5mm
Height	200mm	
Width	65mm	
Net Weight	580g	
Accessory	Penlight	

***MADE IN JAPAN**

MICROSCOPE SET

Biological microscope 40x,100x and 400x, built-in LED illuminator with mirror, dissecting set and prepared slides in a box.



Model No.	M8070
Eyepiece	10x
Objectives	4x, 10x and 40x
Focusing Adjustment	Stage up-and-down movement type
Function	Inclined angle 45 deg.
Stage	90x90mm
Diaphragm	Rotatable disc type (8steps)
Illuminator	LED with mirror(Battery : 3 pcs. AA)
Height / Width	250x105mm
Weight	350g

STEREO MICROSCOPES

Compact size stereo microscope 20x, built-in LED illuminator, in a carrying case.



M8075

M8076

Model No.	M8075 / M8076
Eyepiece	WF10x(Field number : 18)
Objectives	2x
Focusing Adjustment	Body up-and-down movement type
Base	135x98mm
Illuminator	LED (Battery : 2pcs.AA)
Height	230mm
Weight	570g

POCKET MICROSCOPES WITH LIGHT



M783/M784

Model No.	M783	M784
Magnification	40x	100x
Weight	70g	

RADIUS GAUGE

To measure the lens diopter
Refractive index of radius gauge : 1.525



To calculate the radius of curvature

Model No.	591
Measuring Range	0 - 17D
M/M graduation	0.25D
Weight	80g

***MADE IN JAPAN**